

Mosquito Salivary Gland Dissections

Alida Coppi & Photini Sinnis

A. Solutions

▸ For dissection:

1. 70% ethanol

70% EtOH

For 50 mls:
35 ml 100% EtOH
15 ml ddH₂O

* Use high quality 200 proof ethanol; for example Sigma E702-3 for molecular biology

2. 1× PBS pH 7.4

137.0 mM NaCl
2.7 mM KCl
1.4 mM KH₂PO₄
10.0 mM Na₂HPO₄

For 1 L:
8.00 g NaCl
0.20 g KCl
0.24 g KH₂PO₄ (anhydrous)
1.44 g Na₂HPO₄ (anhydrous)
↑ 1 L ddH₂O

* pH should be around 7.4. Autoclave or filter sterilize solution.

* This solution can be made as a 10× stock.

3. *Aedes* saline pH 7.00

150.0 mM NaCl
4.0 mM KCl
3.0 mM CaCl₂
1.8 mM NaHCO₃
0.6 mM MgCl₂
25.0 mM HEPES

For 500 ml:
4.40 g NaCl
0.15 g KCl
0.22 g CaCl₂•2H₂O
0.08 g NaHCO₃
0.06 g MgCl₂•6H₂O
3.00 g HEPES
↑ 500 ml ddH₂O

* pH should be around 7.0. Autoclave or filter sterilize solution.

B. Protocol

1. Obtain female mosquitoes from cages and transfer to 14 ml polypropylene round bottom tubes (Falcon #35-2059).
2. Place tube on ice to anaesthetize the mosquitoes (~ 10 – 30 min.)


3. While mosquitoes are on ice, prepare:
 - a. dissecting microscope with a glass plate
 - b. 2 × 25g needles (one slightly bent) attached to 1 cc syringes
 - c. 2 small Petri dishes (polystyrene 60 × 15 mm Corning #25010):
 - 1 with 70% ethanol; 1 with 1× PBS pH 7.4
 - d. *Aedes* saline (dissection medium depends upon experiment to be performed - can use RPMI or DMEM as well)
4. Place anaesthetized mosquitoes in Petri dish containing 70% ethanol and gently swirl the ethanol to mix (so to completely saturate the mosquitoes), about 20 – 30 sec.

* Do not leave mosquitoes in ethanol for a long period of time as the mosquitoes become soft and difficult to dissect.
5. With curved serrated micro-dissecting forceps, transfer the mosquitoes to Petri dish containing 1× PBS pH 7.4.

* If the dissection will take a while, do not leave mosquitoes in saline for a long period of time as the mosquitoes become soft and difficult to dissect. Instead wet the mosquitoes in batches of no more than 50 at a time.
6. Place a small amount of *Aedes* saline on the glass slide (step 3a) (~ ≤100 μl).
7. With curved serrated micro-dissecting forceps line mosquitoes on the glass plate near the *Aedes* saline.


NOTE: It is easier to dissect if mosquitoes are on their side with their heads upward.

Example:



8. Using 2 syringes with 25g needles, dissect salivary glands (see diagram below):
 - a. With the bent needle (in the left hand), press on the thorax.

- b. Place the other needle (in the right hand) under the head of the mosquito and pull upward.
- c. The salivary glands will be found between the head and thorax.
- d. Cut the material above and below the salivary glands and move the glands to the *Aedes* saline. It is important to keep the salivary glands hydrated.



* If the number of dissections is great, the salivary glands will begin to stick to needles. Therefore, change the needles as is necessary to prevent loss of salivary glands.

9. Using a long Pasteur pipette transfer the medium containing the salivary glands to a tube and keep on ice until dissection is finished or the glands are ready to be used.

* Do not let the glands enter the wide part of the Pasteur pipet because it is difficult to get them out.

* A pulled Pasteur pipette can also be used to transfer the salivary glands. The salivary glands will stick to the sharp end of the pipette.

10. To pellet glands, spin for 5 sec at $16,000 \times g$ (14000 rpm using Eppendorf 5415c) at 4°C . *Aedes* saline can be removed and replaced with fresh saline or other media.

11. Release sporozoites from glands by grinding with a small homogenizer

⇒ If salivary glands are in a large volume of dissection medium:

- grind glands in a Dounce homogenizer

OR

- pellet glands by spinning for 5 sec at $16,000 \times g$ (14000 rpm using Eppendorf 5415c) at 4°C ; dissection medium can be removed for a smaller total volume.

⇒ A greater number of sporozoites can be counted if tube is left on ice a few minutes after homogenization.

For homogenization and counting sporozoites - see next protocol